


DNAJB4 Antibody (F53879)

Catalog No.	Formulation	Size
F53879-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53879-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9UDY4
Applications	Western Blot : 1:1000-2000
Limitations	This DNAJB4 antibody is available for research use only.



Western blot testing of human HeLa cell lysate with DNAJB4 antibody at 1:2000.
Predicted molecular weight: 38 kDa.

Description

DnaJ homolog subfamily B member 4 is likely a chaperone protein.

Application Notes

Titration of the DNAJB4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 216-248 from the human protein was used as the immunogen for the DNAJB4 antibody.

Storage

Aliquot the DNAJB4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.