


Dnajc5 Antibody (R36261)

Catalog No.	Formulation	Size
R36261-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80331
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Dnajc5 antibody is available for research use only.



Western blot testing of human cerebellum lysate with Dnajc5 antibody at 0.1ug/ml.
Predicted molecular weight: 22 kDa.

Description

Additional name(s) for this target protein: DnaJ (Hsp40) homolog, subfamily C, member 5

Application Notes

Optimal dilution of the Dnajc5 antibody should be determined by the researcher.

Immunogen

Amino acids TQLTADSHPSYHTD were used as the immunogen for this Dnajc5 antibody.

Storage

Aliquot and store the Dnajc5 antibody at -20oC.