


## DOCK1 Antibody (R36308)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R36308-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                         |
|----------------------|---------------------------------------------------------|
| Availability         | 1-3 business days                                       |
| Species Reactivity   | Human                                                   |
| Predicted Reactivity | Cow, Dog, Mouse, Pig, Rat                               |
| Format               | Antigen affinity purified                               |
| Clonality            | Polyclonal (goat origin)                                |
| Isotype              | Goat Ig                                                 |
| Purity               | Antigen affinity                                        |
| Gene ID              | 1793                                                    |
| Applications         | Western Blot : 4-8ug/ml                                 |
| Limitations          | This DOCK1 antibody is available for research use only. |



Western blot testing of human A431 lysate with DOCK1 antibody at 4ug/ml. Predicted molecular weight: 215kDa, routinely observed at 180–215kDa.

### Description

Additional name(s) for this target protein: Dedicator of cyto-kinesis 1; DOCK180

### Application Notes

Optimal dilution of the DOCK1 antibody should be determined by the researcher.

### Immunogen

Amino acids TTRKQTSVDSGIVQ were used as the immunogen for this DOCK1 antibody.

## Storage

Aliquot and store the DOCK1 antibody at -20oC.