


DOPA decarboxylase Antibody (R36345)

Catalog No.	Formulation	Size
R36345-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1644
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This DOPA decarboxylase antibody is available for research use only.



Western blot testing of human kidney lysate with DOPA decarboxylase antibody at 0.03ug/ml. Predicted molecular weight: ~54kDa.

Description

Additional name(s) for this target protein: DDC; aromatic L-amino acid decarboxylase (AADC)

Application Notes

Optimal dilution of the DOPA decarboxylase antibody should be determined by the researcher.

Immunogen

Amino acids WEHIKELAADV L were used as the immunogen for this DOPA decarboxylase antibody.

Storage

Aliquot and store the DOPA decarboxylase antibody at -20oC.