


## DPF2 Antibody / UBID4 (R36460)

Catalog No.	Formulation	Size
R36460-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q92785
Gene ID	5977
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This DPF2 antibody is available for research use only.



Western blot testing of human 1) Jurkat and 2) K562 cell lysate with DPF2 antibody at 1ug/ml. Predicted molecular weight: ~44 kDa, routinely observed at ~50 kDa.

## Description

Additional name(s) for this target protein: D4, zinc and double PHD fingers family 2; UBID4; BRG1-associated factor 45D; BRAF45D

## Application Notes

Optimal dilution of the DPF2 antibody should be determined by the researcher.

## Immunogen

Amino acids 245-259 (QPPTPVSQRSEEQKS) were used as the immunogen for this DPF2 antibody.

## **Storage**

Aliquot and store the DPF2 antibody at -20oC.