


Dynactin 3 Antibody (R35494)


Catalog No.	Formulation	Size
R35494-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11258
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Dynactin 3 antibody is available for research use only.



Western blot testing of human kidney lysate with Dynactin 3 antibody at 0.03ug/ml. Expected molecular weight (human): ~22kDa.



Western blot of mouse (A) and rat (B) skeletal muscle lysates with Dynactin 3 antibody at 0.2ug/ml. Expected molecular weight (mouse): ~24kDa.

Description

Additional name(s) for this target protein: DCTN3, DCTN-22, DCTN22; Human: Dynactin 3 (p22); Mouse: Dynactin 3 (p24)

Application Notes

Optimal dilution of the Dynactin 3 antibody should be determined by the researcher.

Immunogen

Amino acids QVALGNISSKRE were used as the immunogen for this Dynactin 3 antibody.

Storage

Aliquot and store the Dynactin 3 antibody at -20oC.