


Dynactin Antibody (R35253)


Catalog No.	Formulation	Size
R35253-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1639
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 2-4ug/ml Immunocytochemistry : suitable ELISA (peptide) LOD : 1:128000
Limitations	This Dynactin antibody is available for research use only.



Western blot testing of human testis lysate with Dynactin antibody at 1ug/ml. Expected molecular weight: 142-150 kDa.



IHC staining of FFPE human cerebellum with Dynactin antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Dynactin1, DCTN1, DCTN

Application Notes

Optimal dilution of the Dynactin antibody should be determined by the researcher.

Immunogen

Amino acids QEQLHQLHSRLIS were used as the immunogen for this Dynactin antibody.

Storage

Aliquot and store the Dynactin antibody at -20oC.