


Dyrk1a Antibody (F47340)

Catalog No.	Formulation	Size
F47340-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47340-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q61214
Applications	Western Blot : 1:1000
Limitations	This Dyrk1a antibody is available for research use only.



Description

Dyrk1a may play a role in a signaling pathway regulating nuclear functions of cell proliferation. Phosphorylates serine, threonine and tyrosine residues in its sequence and in exogenous substrates.

Application Notes

Titration of the Dyrk1a antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 718-745 from the mouse protein was used as the immunogen for this Dyrk1a antibody.

Storage

Aliquot the Dyrk1a antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.