


DYRK1A Antibody (R35024)

Catalog No.	Formulation	Size
R35024-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1859
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This DYRK1A antibody is available for research use only.



Western blot testing of K562 lysate with DYRK1A antibody at 0.5ug/ml. Predicted molecular weight ~85 kDa.

Description

Additional name(s) for this target protein: Dual-specificity tyrosine-(Y)-phosphorylation regulated kinase 1A

Application Notes

Optimal dilution of the DYRK1A antibody should be determined by the researcher.

Immunogen

Amino acids PHSQYSDRRQPN were used as the immunogen for this DYRK1A antibody.

Storage

Aliquot and store the DYRK1A antibody at -20oC.