


DYX1C1 Antibody (R35064)

Catalog No.	Formulation	Size
R35064-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	161582
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This DYX1C1 antibody is available for research use only.



Western blot testing of human brain frontal cortex lysate with DYX1C1 antibody at 0.5ug/ml. Predicted molecular weight: ~48/44/44kDa (isoforms a/b/c).

Description

Additional name(s) for this target protein: Dyslexia susceptibility 1 candidate 1, EKN1

Application Notes

Optimal dilution of the DYX1C1 antibody should be determined by the researcher.

Immunogen

Amino acids PLQVSDYSWQQTKT were used as the immunogen for this DYX1C1 antibody. This sequence is common to

isoforms a, b and c.

Storage

Aliquot and store the DYX1C1 antibody at -20oC.