


## E2F3 Antibody (RQ6616)

| Catalog No. | Formulation                                           | Size   |
|-------------|-------------------------------------------------------|--------|
| RQ6616      | 0.5mg/ml if reconstituted with 0.2ml sterile DI water | 100 ug |

**Bulk quote request**

|                           |                                                        |
|---------------------------|--------------------------------------------------------|
| <b>Availability</b>       | 1-3 business days                                      |
| <b>Species Reactivity</b> | Human, Mouse, Rat                                      |
| <b>Format</b>             | Antigen affinity purified                              |
| <b>Clonality</b>          | Polyclonal (rabbit origin)                             |
| <b>Isotype</b>            | Rabbit IgG                                             |
| <b>Purity</b>             | Antigen affinity purified                              |
| <b>Buffer</b>             | Lyophilized from 1X PBS with 2% Trehalose              |
| <b>UniProt</b>            | O00716                                                 |
| <b>Applications</b>       | Western Blot : 1-2ug/ml                                |
| <b>Limitations</b>        | This E2F3 antibody is available for research use only. |



Western blot testing of 1) human 293T, 2) human SK-O-V3, 3) human K562, 4) human HEL, 5) rat brain and 6) mouse brain tissue lysate with E2F3 antibody. Predicted molecular weight: ~49 kDa but can also be observed from 56~66 kDa.

## Description

Transcription factor E2F3, also known as KIAA0075, is a protein that in humans is encoded by the E2F3 gene. The protein encoded by this gene is a member of the E2F family of transcription factors. By fluorescence in situ hybridization, E2F3 gene is mapped to 6q22.3. The induction of specific E2F activities is an essential component in the MYC pathways that control cell proliferation and cell fate decisions. This gene encodes a member of a small family of transcription factors that function through binding of DP interaction partner proteins. The encoded protein recognizes a specific sequence motif in DNA and interacts directly with the retinoblastoma protein (pRB) to regulate the expression of genes involved in the cell cycle.

## Application Notes

Optimal dilution of the E2F3 antibody should be determined by the researcher.

## Immunogen

Amino acids KTNNQDHNGNIPKPASKDLA from the human protein were used as the immunogen for the E2F3 antibody.

## Storage

After reconstitution, the E2F3 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.