


E2F4 Antibody (R35037)

Catalog No.	Formulation	Size
R35037-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1874
Localization	Cytoplasmic, nuclear
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This E2F4 antibody is available for research use only.



Western blot testing of A431 [A] and HeLa [B] lysates with E2F4 antibody at 0.3ug/ml. Expected molecular weight ~44 kDa (unmodified) and 60-65 kDa (phosphorylated).

Description

Additional name(s) for this target protein: Transcription factor E2F4, E2F transcription factor 4

Application Notes

Optimal dilution of the E2F4 antibody should be determined by the researcher.

Immunogen

Amino acids PKELSEIFDPTR were used as the immunogen for this E2F4 antibody.

Storage

Aliquot and store the E2F4 antibody at -20oC.