


EEA1 Antibody (R35258)

Catalog No.	Formulation	Size
R35258-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8411
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This EEA1 antibody is available for research use only.



Western blot testing of HEK293 lysate with EEA1 antibody at 0.3ug/ml. Predicted molecular weight ~162kDa.

Description

Additional name(s) for this target protein: Early endosome antigen 1, MST105, MSTP105

Application Notes

Optimal dilution of the EEA1 antibody should be determined by the researcher.

Immunogen

Amino acids QETKIQHEELNNRIQ were used as the immunogen for this EEA1 antibody.

Storage

Aliquot and store the EEA1 antibody at -20oC.