


EGFL7 Antibody (R34184)

Catalog No.	Formulation	Size
R34184-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51162
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This EGFL7 antibody is available for research use only.



Western blot testing of human heart [A], kidney [B] and lung [C] lysates with EGFL7 antibody at 1ug/ml. Predicted molecular weight: ~29/27kDa (pro/mature form).

Description

Additional name(s) for this target protein: EGF-like-domain, multiple 7

Application Notes

Optimal dilution of the EGFL7 antibody should be determined by the researcher.

Immunogen

Amino acids HGDPVSESFVQR were used as the immunogen for this EGFL7 antibody.

Storage

Aliquot and store the EGFL7 antibody at -20oC.