


Elongator complex protein 3 Antibody / ELP3 / KAT9 [clone 8B4] (RQ6406)

Catalog No.	Formulation	Size
RQ6406	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2a
Clone Name	8B4
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q9H9T3
Localization	Cytoplasmic, nuclear
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Elongator complex protein 3 antibody is available for research use only.



Western blot testing of human 1) K562, 2) U-2 OS, 3) HepG2, 4) U-87 MG, 5) rat lung, 6) rat PC-12, 7) mouse lung and 8) mouse NIH 3T3 cell lysate with Elongator complex protein 3 antibody. Predicted molecular weight ~62 kDa.

Description

Elongator complex protein 3, also named KAT9, is a protein that in humans is encoded by the ELP3 gene. It is mapped to 8p21.1. ELP3 is the catalytic histone acetyltransferase subunit of the RNA polymerase II elongator complex, which is a component of the RNA polymerase II (Pol II) holoenzyme and is involved in transcriptional elongation. ELP3 supports the migration and branching of projection neurons through acetylation of alpha-tubulin in the developing cerebral cortex. In mammals, ELP3 is important for paternal DNA demethylation after fertilization. ELP3 is potentially involved in cellular redox

homeostasis by mediating the acetylation of glucose-6-phosphate dehydrogenase. Besides, ELP3 may play a role in chromatin remodeling and is involved in acetylation of histones H3 and probably H4.

Application Notes

Optimal dilution of the Elongator complex protein 3 antibody should be determined by the researcher.

Immunogen

Amino acids RDPTKFQHQQFGMLLMEEAERIAREEH were used as the immunogen for the Elongator complex protein 3 antibody.

Storage

After reconstitution, the Elongator complex protein 3 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.