


EpCAM Antibody PE Conjugate [clone EPM17-1] (V7010PE)

Catalog No.	Formulation	Size
V7010PE-100T	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 Tests

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	PE Conjugate
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	EPM17-1
Purity	Protein G affinity chromatography
UniProt	P16422
Localization	Cell surface
Applications	Flow Cytometry : 5ul/test/10e6 cells in 0.1ml Immunofluorescence : 1:50-1:100
Limitations	This EpCAM antibody is available for research use only.



Flow cytometry testing of human MCF-7 cells using PE conjugated EpCAM antibody (red) and isotype control (green).

Description

EpCAM may act as a physical homophilic interaction molecule between intestinal epithelial cells (IECs) and intraepithelial lymphocytes (IELs) at the mucosal epithelium for providing immunological barrier as a first line of defense against mucosal infection. [UniProt]

Application Notes

Titering of the EpCAM antibody may be required for optimal performance.

Immunogen

Recombinant human protein was used as the immunogen for the EpCAM antibody.

Storage

Store the EpCAM antibody at 2-8°C. Conjugate is light sensitive, store in the dark.