


## Ephrin B2 Antibody (R33300)

Catalog No.	Formulation	Size
R33300-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Predicted Reactivity</b>	Cow, Dog, Pig
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	1948
<b>Applications</b>	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
<b>Limitations</b>	This Ephrin B2 antibody is available for research use only.



Western blot testing of rat lung lysate with Ephrin B2 antibody at 0.1ug/ml. Predicted molecular weight ~37kDa.

### Description

Additional name(s) for this target protein: EFNB2

### Application Notes

Optimal dilution of the Ephrin B2 antibody should be determined by the researcher.

### Immunogen

Amino acids CPKVDSKTVGQYE were used as the immunogen for this Ephrin B2 antibody.

## Storage

Aliquot and store the Ephrin B2 antibody at -20oC.