


ERAB Antibody / HSD17B10 (R33254)

Catalog No.	Formulation	Size
R33254-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3028
Applications	Western Blot : 0.3-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ERAB antibody is available for research use only.



Western blot testing of human brain lysate with ERAB antibody at 0.3ug/ml. Predicted molecular weight: ~27kDa.

Description

Additional name(s) for this target protein: Hydroxysteroid (17-beta) dehydrogenase 10; HSD17B10

Application Notes

Optimal dilution of the ERAB antibody should be determined by the researcher.

Immunogen

Amino acids CIRLDGAIRMQP were used as the immunogen for this ERAB antibody.

Storage

Aliquot and store the ERAB antibody at -20oC.