


ErbB3 Antibody (F53392)

Catalog No.	Formulation	Size
F53392-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53392-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q61526
Applications	Western Blot : 1:2000
Limitations	This ErbB3 antibody is available for research use only.



Western blot testing of ErbB3 antibody at 1:2000 dilution + mouse lung lysate; Predicted band size : 148 kDa.

Description

Binds and is activated by neuregulins and NTAK. May also be activated by CSPG5.

Application Notes

Titration of the ErbB3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This mouse ErbB3 antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 1076-1110 amino acids from the region of mouse ErbB3.

Storage

Aliquot the Erbb3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.