


EVL Antibody (R36026)

Catalog No.	Formulation	Size
R36026-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51466
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This EVL antibody is available for research use only.



Western blot testing of Jurkat lysate with EVL antibody at 1ug/ml. Predicted molecular weight: ~45kDa.

Description

Additional name(s) for this target protein: Enah/Vasp-like

Application Notes

Optimal dilution of the EVL antibody should be determined by the researcher.

Immunogen

Amino acids SQPPNSSEAGRKP were used as the immunogen for this EVL antibody.

Storage

Aliquot and store the EVL antibody at -20oC.