


FABP2 Antibody Biotin Conjugate (R33789BTN)

Catalog No.	Formulation	Size
R33789BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2169
Localization	Cytoplasmic
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:64000
Limitations	This FABP2 antibody is available for research use only.



Western blot of mouse duodenum lysate with biotinylated FABP2 antibody at 0.1ug/ml.
Predicted molecular weight: ~15kDa.

Description

Additional name(s) for this target protein: Fatty acid binding protein 2

Application Notes

Optimal dilution of the FABP2 antibody should be determined by the researcher.

Immunogen

Amino acids EGVEAKRIFKKD were used as the immunogen for this FABP2 antibody.

Storage

Aliquot and store the FABP2 antibody at -20oC.