


FABP2 Antibody (R33789)


Catalog No.	Formulation	Size
R33789-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2169
Localization	Cytoplasmic
Applications	Western Blot : 0.001-0.1ug/ml Immunohistochemistry (FFPE) : 2-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This FABP2 antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing FABP2 with FABP2 antibody [mock transfection in lane 1]. Predicted molecular weight: ~15 kDa.



IHC staining of FFPE human small intestine with FABP2 antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.



Western blot of human duodenum lysate with FABP2 antibody at 0.001ug/ml. Predicted molecular weight: ~15 kDa.

Description

Additional name(s) for this target protein: Fatty acid binding protein 2

Application Notes

Optimal dilution of the FABP2 antibody should be determined by the researcher.

Immunogen

Amino acids EGVEAKRIFKKD were used as the immunogen for this FABP2 antibody.

Storage

Aliquot and store the FABP2 antibody at -20oC.