


## FABP3 Antibody (cardiac) (R35172)

Catalog No.	Formulation	Size
R35172-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2170
Localization	Cytoplasmic, membrane
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FABP3 antibody is available for research use only.



Western blot testing of human [A], mouse [B] and rat [C] heart lysates with FABP3 antibody at 0.01ug/ml. Predicted molecular weight: ~15 kDa.

### Description

Additional name(s) for this target protein: Fatty acid binding protein 3, heart-type FABP

### Application Notes

Optimal dilution of the FABP3 antibody should be determined by the researcher.

## Immunogen

Amino acids PTIIIEKNGDILTLK were used as the immunogen for this FABP3 antibody.

## Storage

Aliquot and store the FABP3 antibody at -20oC.