


Factor XIIIa Antibody (R34270)

Catalog No.	Formulation	Size
R34270-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2162
Applications	Western Blot : 0.05-0.2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Factor XIIIa antibody is available for research use only.



Western blot testing of human placenta lysate with Factor XIIIa antibody at 0.1ug/ml.
Predicted molecular weight: ~83 kDa.

Description

Additional name(s) for this target protein: Coagulation factor XIII A1 polypeptide; F13A1

Application Notes

Optimal dilution of the Factor XIIIa antibody should be determined by the researcher.

Immunogen

Amino acids HRKLIASMSSDSL RH were used as the immunogen for this Factor XIIIa antibody.

Storage

Aliquot and store the Factor XIIIa antibody at -20oC.