


FAM129B Antibody / MINERVA / NIBAN2 (RQ6569)


Catalog No.	Formulation	Size
RQ6569	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q96TA1
Localization	Cytoplasmic
Applications	Western Blot : 1-2ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml Flow Cytometry : 1-3ug/million cells Direct ELISA : 0.1-0.5ug/ml
Limitations	This FAM129B antibody is available for research use only.




IHC staining of FFPE human lung cancer tissue with FAM129B antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE human placental tissue with FAM129B antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of 1) human HeLa, 2) human RT4, 3) rat stomach, 4) rat C6, 5) mouse stomach and 6) mouse NIH 3T3 cell lysate with FAM129B antibody. Predicted molecular weight ~84 kDa.



Flow cytometry testing of human HeLa cells with FAM129B antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= FAM129B antibody.

Description

FAM129B/Niban-like protein 1 (family with sequence similarity 129, member B) belongs to a poorly characterized family of Niban proteins that also includes FAM129A/Niban and FAM129C/Niban-like protein 2. FAM129A/Niban is implicated in the ER stress response and is upregulated at the protein level in thyroid carcinoma. FAM129C/Niban-like protein 2 is preferentially expressed in B-cells and is one of five biomarkers upregulated in whole blood from patients with colorectal carcinoma. FAM129B is broadly expressed and has been shown to be a downstream target of B-Raf in melanoma cells. Though FAM129B does not appear to regulate cell growth and division, phosphorylation of FAM129B by B-Raf is essential for the invasive potential of melanoma and non-melanoma cell lines. Deletion of FAM129B in melanoma cells significantly impairs B-Raf/MEK/Erk-dependent invasion into the extracellular matrix.

Application Notes

Optimal dilution of the FAM129B antibody should be determined by the researcher.

Immunogen

An E. coli-derived human protein (amino acids K23-K561) was used as the immunogen for the FAM129B antibody.

Storage

After reconstitution, the FAM129B antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.

