


Fascin-3 Antibody (F49847)


Catalog No.	Formulation	Size
F49847-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F49847-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q9NQT6
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:50-1:100 Flow Cytometry : 1:10-1:50
Limitations	This Fascin-3 antibody is available for research use only.



Fascin-3 antibody flow cytometric analysis of MDA-MB231 cells (green) compared to a [negative control](http://search_result.php?search_txt=n1001) (blue). FITC-conjugated goat-anti-rabbit secondary Ab was used for the analysis.



Fascin-3 antibody IHC analysis in formalin fixed and paraffin embedded human testis tissue.



Western blot analysis of Fascin-3 antibody and 1) MDA-MB231, 2) NCI-H460 cell line and 3) mouse testis tissue lysate.

Description

FSCN3 acts as an actin bundling protein.

Application Notes

Titration of the Fascin-3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-30 from the human protein was used as the immunogen for this Fascin-3 antibody.

Storage

Aliquot the Fascin-3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.