

FGF1 Antibody (R30107)


Catalog No.	Formulation	Size
R30107	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 0.025% sodium azide
UniProt	P05230
Localization	Cytoplasmic
Applications	Western Blot (recombinant Protein) : 0.5-1ug/ml IHC (FFPE) : 0.5-1ug/ml ELISA : 0.1-0.5ug/ml (human protein tested)
Limitations	This FGF1 antibody is available for research use only.



Western blot testing of FGF1 antibody and Lane 1: 10ng; 2: 5ng; 3: 2.5ng; 4: SMMC-7721 cell lysate. Expected molecular weight ~17 kDa.



IHC-P: FGF1 antibody testing of human breast cancer tissue

Description

Fibroblast growth factor 1 (acidic), also known as [FGF1](#)/ECGF/HBGF1, is a human gene which is mapped to 5q31. Human FGF1 shares 96% amino acid sequence homology with both rat and mouse. The protein encoded by this gene is a member of the fibroblast growth factor family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This protein functions as a modifier of endothelial cell migration and proliferation, as well as an angiogenic factor. It acts as a mitogen for a variety of mesoderm- and neuroectoderm-derived cells in vitro, thus is thought to be involved in organogenesis. Additionally, Acidic fibroblast growth factor is derived from beta-endothelial cell growth factor (ECGFB) by posttranslational processing. Alpha-ECGF is also derived from ECGFB in the same manner.

Application Notes

The stated application concentrations are suggested starting points. Titration of the FGF1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Human partial recombinant protein (AA 16-155) was used as the immunogen for this FGF1 antibody.

Storage

After reconstitution, the FGF1 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.