


FGF21 Antibody (R33314)

Catalog No.	Formulation	Size
R33314-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26291
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This FGF21 antibody is available for research use only.



Western blot testing of human liver lysate with FGF21 antibody at 1ug/ml. Predicted molecular weight ~21kDa.

Description

Additional name(s) for this target protein: Fibroblast growth factor 21

Application Notes

Optimal dilution of the FGF21 antibody should be determined by the researcher.

Immunogen

Amino acids CQRPDGALYGSLH were used as the immunogen for this FGF21 antibody.

Storage

Aliquot and store the FGF21 antibody at -20oC.