


Fgfr1 Antibody (F44216)


Catalog No.	Formulation	Size
F44216-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44216-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P16092
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This Fgfr1 antibody is available for research use only.



Fgfr1 antibody western blot analysis in NCI-H460 lysate. Predicted molecular weight: 75-160 kDa depending on glycosylation level.



Fgfr1 antibody immunohistochemistry analysis in formalin fixed and paraffin embedded mouse adrenal glands.

Description

Receptor for basic fibroblast growth factor. A shorter form of the receptor could be a receptor for FGF1 (aFGF). Receptor for FGF23 in the presence of KL.

Application Notes

Titration of the Fgfr1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 756-784 from the mouse protein was used as the immunogen for this Fgfr1 antibody.

Storage

Aliquot the Fgfr1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.