


FGFR2 Antibody (R34123)

Catalog No.	Formulation	Size
R34123-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2263
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This FGFR2 antibody is available for research use only.



Western blot testing of A549 lysate with FGFR2 antibody at 0.3ug/ml. Predicted molecular weight of multiple isoforms: 80-120 kDa. The observed size may be larger due to glycosylation.

Description

Additional name(s) for this target protein: Fibroblast growth factor receptor 2

Application Notes

Optimal dilution of the FGFR2 antibody should be determined by the researcher.

Immunogen

Amino acids GREKEITASPDY were used as the immunogen for this FGFR2 antibody. The amino acid sequence used as

immunogen is from the extracellular portion of the protein.

Storage

Aliquot and store the FGFR2 antibody at -20oC.