


FGFR3 Antibody (F44217)


Catalog No.	Formulation	Size
F44217-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44217-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q61851
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This FGFR3 antibody is available for research use only.



FGFR3 antibody immunohistochemistry analysis in formalin fixed and paraffin embedded mouse brain tissue.



FGFR3 antibody western blot analysis in 293 lysate. Predicted molecular weight: 87-135 kDa depending on glycosylation level.

Description

FGFR3 is a receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.

Application Notes

Titration of the FGFR3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 468-495 from the mouse protein was used as the immunogen for this FGFR3 antibody.

Storage

Aliquot the FGFR3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.