


FLI1 Antibody (F49648)


Catalog No.	Formulation	Size
F49648-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F49648-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q01543
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:50-1:100
Limitations	This FLI1 antibody is available for research use only.



Western blot analysis of FLI1 antibody and 293 lysate.



FLI1 antibody IHC analysis in formalin fixed and paraffin embedded human tonsil tissue.

Description

FLI1 is a sequence-specific transcriptional activator. It recognizes the DNA sequence 5'-C[CA]GGAAGT-3'.

Application Notes

Titration of the FLI1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 52-80 from the human protein was used as the immunogen for this FLI1 antibody.

Storage

Aliquot the FLI1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.