


Flotillin 1 Antibody [clone BHG-6] (RQ4792)

Catalog No.	Formulation	Size
RQ4792	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	BHG-6
Purity	Affinity purified
UniProt	O75955
Applications	Western Blot : 1:500-1:2000
Limitations	This Flotillin 1 antibody is available for research use only.



Western blot testing of human HeLa cell lysate with Flotillin 1 antibody. Predicted molecular weight ~49 kDa.

Description

Flotillin 1 may act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles. [UniProt]

Application Notes

Optimal dilution of the Flotillin 1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human Flotillin 1 / FLOT1 was used as the immunogen for the Flotillin 1 antibody.

Storage

Store the Flotillin 1 antibody at -20oC.