


Flotillin 1 Antibody (R35932)


Catalog No.	Formulation	Size
R35932-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10211
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Flotillin 1 antibody is available for research use only.



Western blot testing of H460 lysate with Flotillin 1 antibody at 0.3ug/ml.
Expected/observed molecular weight ~49kDa.



Western blot of HEK293 lysate overexpressing FLOT1 probed with Flotillin 1 antibody (mock transfection in lane 1). Expected/observed molecular weight ~49kDa.



IHC staining of FFPE human kidney with Flotillin antibody at 3.8ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: FLOT1

Application Notes

Optimal dilution of the Flotillin 1 antibody should be determined by the researcher.

Immunogen

Amino acids SISQVNHKPLRTA were used as the immunogen for this Flotillin 1 antibody.

Storage

Aliquot and store the Flotillin 1 antibody at -20oC.