


## FLVCR Antibody (R34680)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R34680-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                               |
|--------------------|---------------------------------------------------------------|
| Availability       | 1-3 business days                                             |
| Species Reactivity | Human                                                         |
| Format             | Antigen affinity purified                                     |
| Clonality          | Polyclonal (goat origin)                                      |
| Isotype            | Goat Ig                                                       |
| Purity             | Antigen affinity                                              |
| Gene ID            | 28982                                                         |
| Applications       | Western Blot : 0.1-0.3ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations        | This FLVCR antibody is available for research use only.       |



Western blot testing of human bone marrow lysate with FLVCR antibody at 0.1ug/ml. The expected ~55kDa band and the additional ~26kDa band are blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: Feline leukemia virus subgroup C cellular receptor; FLVCR1

## Application Notes

Optimal dilution of the FLVCR antibody should be determined by the researcher.

## Immunogen

Amino acids KTVMLSKQSESAI were used as the immunogen for this FLVCR antibody.

## Storage

Aliquot and store the FLVCR antibody at -20oC.