


FOXA2 Antibody (R34150)


Catalog No.	Formulation	Size
R34150-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3170
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 4-8ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This FOXA2 antibody is available for research use only.



Western blot testing of human liver lysate with FOXA2 antibody at 1ug/ml



IHC staining of FFPE human prostate with FOXA2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.



Description

Additional name(s) for this target protein: Forkhead box A2; HNF3B

Application Notes

Optimal dilution of the FOXA2 antibody should be determined by the researcher.

Immunogen

Amino acids GVYSRPIMNSS were used as the immunogen for this FOXA2 antibody.

Storage

Aliquot and store the FOXA2 antibody at -20oC.