


## FOXC1 Antibody (R33387)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33387-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                          |
|----------------------|----------------------------------------------------------|
| Availability         | 1-3 business days                                        |
| Species Reactivity   | Human                                                    |
| Predicted Reactivity | Mouse, Pig                                               |
| Format               | Antigen affinity purified                                |
| Clonality            | Polyclonal (goat origin)                                 |
| Isotype              | Goat Ig                                                  |
| Purity               | Antigen affinity                                         |
| Gene ID              | 2296                                                     |
| Applications         | Western Blot : 2-4ug/ml<br>ELISA (peptide) LOD : 1:32000 |
| Limitations          | This FOXC1 antibody is available for research use only.  |



Western blot testing of HEK293 lysate with FOXC1 antibody at 2ug/ml. Predicted molecular weight: ~57 kDa, routinely observed at up to 75 kDa.

## Description

Additional name(s) for this target protein: Forkhead box C1

## Application Notes

Optimal dilution of the FOXC1 antibody should be determined by the researcher.

## Immunogen

Amino acids DAVKDKEEKDRLH were used as the immunogen for this FOXC1 antibody.

## Storage

Aliquot and store the FOXC1 antibody at -20oC.