


FOXC1 Antibody (R35797)


Catalog No.	Formulation	Size
R35797-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2296
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 4-6ug/ml ChIP : suitable ELISA (peptide) LOD : 1:4000
Limitations	This FOXC1 antibody is available for research use only.



Western blot of HEK293 lysate over expressing human FOXC1-FLAG tested with FOXC1 antibody [1ug/ml] in Lane A anti-FLAG (1/3000) in lane C. Mock-transfected HEK293 tested with FOXC1 antibody [1ug/ml] in Lane B. Predicted molecular weight: ~57 kDa, routinely observed at up to 75 kDa.



IHC staining of FFPE human kidney with FOXC1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Forkhead box C1

Application Notes

Optimal dilution of the FOXC1 antibody should be determined by the researcher.

Immunogen

Amino acids RTSGAFVYDCSKF were used as the immunogen for this FOXC1 antibody.

Storage

Aliquot and store the FOXC1 antibody at -20oC.