


FOXG1 Antibody (R34258)

Catalog No.	Formulation	Size
R34258-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2290
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This FOXG1 antibody is available for research use only.



Western blot testing of human brain lysate with FOXG1 antibody at 0.3ug/ml. Predicted molecular weight ~50/52 kDa (FOXG1A/FOXG1B). Both observed bands are blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Forkhead box G1, BF1, BF2, FHLK3, FKH2

Application Notes

Optimal dilution of the FOXG1 antibody should be determined by the researcher.

Immunogen

Amino acids HQNQGSSSNPLIH were used as the immunogen for this FOXG1 antibody.

Storage

Aliquot and store the FOXG1 antibody at -20oC.