


FOXL2 Antibody (R33635)

Catalog No.	Formulation	Size
R33635-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	668
Localization	Nuclear
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This FOXL2 antibody is available for research use only.



Western blot testing of human ovary lysate with FOXL2 antibody at 1ug/ml. Predicted molecular weight ~38 kDa, but can be observed at 45-50 kDa.

Description

Additional name(s) for this target protein: Forkhead box L2

Application Notes

Optimal dilution of the FOXL2 antibody should be determined by the researcher.

Immunogen

Amino acids DSKTGALHSRLDL were used as the immunogen for this FOXL2 antibody.

Storage

Aliquot and store the FOXL2 antibody at -20oC.

References (1)