


FOXP1 Antibody [clone AOGO-6] (RQ7562)


Catalog No.	Formulation	Size
RQ7562	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AOGO-6
UniProt	Q9H334
Localization	Nuclear
Applications	Immunofluorescence : 1:25-1:50 Western Blot : 1:250-1:500
Limitations	This FOXP1 antibody is available for research use only.



Immunofluorescent staining of FFPE human HeLa cells with FOXP1 antibody (red) and Beta Tubulin antibody (green). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of 1) human Jurkat, 2) human A549, 3) human MCF7, 4) rat brain, 5) mouse spleen and 6) mouse brain tissue cell lysate with FOXP1 antibody at 1:500 dilution. Expected molecular weight: 65-77 kDa (multiple isoforms).

Description

FOXP1 may act as a tumor suppressor as it is lost in several tumor types.

Application Notes

The stated application concentrations are suggested starting points. Titration of the FOXP1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide specific to human FOXP1 was used as the immunogen for the FOXP1 antibody.

Storage

Store the FOXP1 antibody at -20oC.