


FOXQ1 Antibody (R34512)

Catalog No.	Formulation	Size
R34512-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	94234
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FOXQ1 antibody is available for research use only.



Western blot testing of H460 lysate with FOXQ1 antibody at 1ug/ml. Predicted molecular weight: ~46kDa.

Description

Additional name(s) for this target protein: Forkhead box Q1, FKHL1

Application Notes

Optimal dilution of the FOXQ1 antibody should be determined by the researcher.

Immunogen

Amino acids KLEVFVPRAAHGD were used as the immunogen for this FOXQ1 antibody.

Storage

Aliquot and store the FOXQ1 antibody at -20oC.