


Fractalkine Antibody (R34679)

Catalog No.	Formulation	Size
R34679-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6376
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Fractalkine antibody is available for research use only.



Western blot testing of human colorectal cancer lysate with Fractalkine antibody at 0.3ug/ml. Predicted molecular weight: ~42kDa.

Description

Additional name(s) for this target protein: Chemokine (C-X3-C motif) ligand 1, ABCD-3, C3Xkine; CX3CL1

Application Notes

Optimal dilution of the Fractalkine antibody should be determined by the researcher.

Immunogen

Amino acids KTSEAPSTQDPST were used as the immunogen for this Fractalkine antibody.

Storage

Aliquot and store the Fractalkine antibody at -20oC.