

## FTL Antibody / Ferritin light chain Biotin Conjugate (R34085BTN)

Catalog No.	Formulation	Size
R34085BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2512
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:128000
Limitations	This FTL antibody is available for research use only.



Western blot testing of rat brain lysate with biotinylated FTL antibody at 1ug/ml.  
Predicted molecular weight: ~20kDa.

### Description

Additional name(s) for this target protein: Ferritin L, Ferritin light chain

### Application Notes

Optimal dilution of the FTL antibody should be determined by the researcher.

### Immunogen

Amino acids GEYLFERLTKHD were used as the immunogen for this FTL antibody.

## **Storage**

Aliquot and store the FTL antibody at -20oC.