


FTO Antibody (R35410)

Catalog No.	Formulation	Size
R35410-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79068
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FTO antibody is available for research use only.



Western blot testing of human pancreas lysate with FTO antibody at 0.03ug/ml.
Predicted molecular weight ~58kDa.

Description

Additional name(s) for this target protein: Fat mass and obesity associated

Application Notes

Optimal dilution of the FTO antibody should be determined by the researcher.

Immunogen

Amino acids QQKPD CRPYWEKDD were used as the immunogen for this FTO antibody.

Storage

Aliquot and store the FTO antibody at -20oC.