


FZD9 Antibody (R36199)

Catalog No.	Formulation	Size
R36199-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8326
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FZD9 antibody is available for research use only.



IHC testing of FFPE human testis with FZD9 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Frizzled 9; Frizzled homolog 9

Application Notes

Optimal dilution of the FZD9 antibody should be determined by the researcher.

Immunogen

Amino acids TKTDPSLENPTH were used as the immunogen for this FZD9 antibody.

Storage

Aliquot and store the FZD9 antibody at -20oC.