


## GAB2 Antibody (R35496)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R35496-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                              |
|----------------------|--------------------------------------------------------------|
| Availability         | 1-3 business days                                            |
| Species Reactivity   | Human                                                        |
| Predicted Reactivity | Dog, Mouse, Rat                                              |
| Format               | Antigen affinity purified                                    |
| Clonality            | Polyclonal (goat origin)                                     |
| Isotype              | Goat Ig                                                      |
| Purity               | Antigen affinity                                             |
| Gene ID              | 9846                                                         |
| Applications         | Western Blot : 0.5-1.5ug/ml<br>ELISA (peptide) LOD : 1:64000 |
| Limitations          | This GAB2 antibody is available for research use only.       |



Western blot testing of K562 cell lysate with GAB2 antibody at 0.5ug/ml. Predicted molecular weight: ~75kDa.

## Description

Additional name(s) for this target protein: GRB2-associated binding protein 2

## Application Notes

Optimal dilution of the GAB2 antibody should be determined by the researcher.

## Immunogen

Amino acids QVDKEKTQALQN were used as the immunogen for this GAB2 antibody.

## Storage

Aliquot and store the GAB2 antibody at -20oC.