


GABA A Receptor gamma 2 Antibody / GABRG2 (R36452)

Catalog No.	Formulation	Size
R36452-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P18507
Gene ID	2566
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This GABA A Receptor gamma 2 antibody is available for research use only.



Western blot testing of 1) human cerebellum, 2) mouse brain, 3) rat brain and 4) pig brain lysate with GABA A Receptor gamma 2 antibody at 0.1ug/ml. Predicted molecular weight: ~54 kDa.

Description

Additional name(s) for this target protein: GABRG2; Gamma-aminobutyric acid receptor subunit gamma-2

Application Notes

Optimal dilution of the GABA A Receptor gamma 2 antibody should be determined by the researcher.

Immunogen

Amino acids 384-396 (KAPTIDIRPRSAT) were used as the immunogen for this GABA A Receptor gamma 2 antibody.

Storage

Aliquot and store the GABA A Receptor gamma 2 antibody at -20oC.