


## GABRG2 Antibody (R34865)

Catalog No.	Formulation	Size
R34865-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2566
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This GABRG2 antibody is available for research use only.



Western blot testing of mouse brain lysate with GABRG2 antibody at 2ug/ml. Predicted molecular weight: ~55kDa.

## Description

Additional name(s) for this target protein: Gamma-aminobutyric acid (GABA) A receptor, gamma 2

## Application Notes

Optimal dilution of the GABRG2 antibody should be determined by the researcher.

## Immunogen

Amino acids NNATHLQERDE were used as the immunogen for this GABRG2 antibody.

## Storage

Aliquot and store the GABRG2 antibody at -20oC.