


GAL1 Antibody / Galectin 1 [clone 28L34] (RQ8282)

Catalog No.	Formulation	Size
RQ8282	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant RABBIT MONOCLONAL

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28L34
Purity	Affinity chromatography
UniProt	P09382
Applications	Western Blot : 1:500-1:2000
Limitations	This GAL1 antibody is available for research use only.



Western blot testing of human HEK293 cell lysate with GAL1 antibody. Predicted molecular weight ~14 kDa.

Description

Lectin that binds beta-galactoside and a wide array of complex carbohydrates. Plays a role in regulating apoptosis, cell proliferation and cell differentiation. [UniProt]

Application Notes

Optimal dilution of the GAL1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide derived from human Galectin 1 was used as the immunogen for the GAL1 antibody.

Storage

Store the GAL1 antibody at -20oC.